

Product Datasheet

GX5533 - Fast Panoptic Staining Kit for Haematology (RUO)

Product Details

Product Name	Fast Panoptic Staining Kit for Haematology (RUO)
Glenthams Code	GX5533
CAS Number	-
EINECS	-

Structure

Molecular Weight	: -
Molecular Formula	: -



Storage

Recommended storage temperature: +20°C.

Hazards and Transport

CLP Classification	Acute Tox. 3, Acute Tox. 3, Flam. Liq. 2, STOT SE 1, Acute Tox. 3
Signal Word	Danger
Hazard Codes	H311, H331, H225, H370, H301

Precautionary Codes

P280, P303+P361+P353,
P210, P304+P340,
P403+P233, P403+P235,
P261, P264

Pictograms



UN Number, Class, PG (RID)

UN1230 - 3 (6.1) - PG II

Glenthams Product Specification

Suitability for microscopy	: To pass test
Nuclei	: Purple to blue
Lymphocyte cytoplasm	: Blue
Neutrophilic granules	: Violet
Basophilic granules	: Dark violet
Monocyte cytoplasm	: Pink to grey
Eosinophilic granules	: Red to violet
Erythrocytes	: Pink
Platelets	: Purple
Version	: v1.0

About Fast Panoptic Staining Kit for Haematology (RUO)

The Fast Panoptic Staining Kit enables quick and efficient staining of haematological samples for microscopic analysis using a modified Romanowsky-type staining method. This kit is for Research Use Only (RUO).

Two variants of the kit are available: 3 x 100ml for up to 250 applications, and 3 x 500ml for up to 1250 applications.

This document was generated electronically and is therefore valid without signature. © Glenthams Life Sciences Ltd, 2025