Glentham LIFE SCIENCES

Glentham Life Sciences Ltd Unit 5 Leafield Way Corsham SN13 9SW United Kingdom

- t: +44 (0) 1225 667 798
- f: +44 (0) 2033 978 909
- e: info@glentham.com

w: www.glentham.com

Product Datasheet

GE1836 - Laemmli Sample Buffer 4X (without 2-Mercaptoethanol)

Product Details

Product Name	Laemmli Sample Buffer 4X (without 2-Mercaptoethanol)
Glentham Code	GE1836
CAS Number	-
Related Categories	Biochemicals, Reagents for Gel Electrophoresis of DNA/RNA, Reagents for Gel Electrophoresis of Proteins, Buffer Solutions

About Laemmli Sample Buffer 4X (without 2-Mercaptoethanol)

Laemmli sample buffer is used in protein sample preparation for the Laemmli system, the most commonly used SDS-PAGE method. GE1836 is a 4X strength solution without 2-Mercaptoethanol.

Structure

Molecular Weight	:	-
Molecular Formula	:	-

Storage

Recommended storage temperature: +4°C.

Hazards and Transport

Not classified as hazardous under CLP. Not classified as dangerous for transport.

Glentham Product Specification

Physical Description	:	Clear, blue to very dark blue liquid
pН	:	~6.8
Composition	:	TRIS: 0.2M, SDS: 0.277M, Glycerol: 40%, Bromophenol Blue: 3.73mM
Version	:	v1.0

This document was generated electronically and is therefore valid without signature. Glentham Life Sciences Ltd, 2025